

### Peptide Amphiphile Micelles Aid in Diagnostics of Atherosclerosis

Iris Hsu, Deborah Chin, Eun Ji Chung
2hsuir2012@gmail.com
Diamond Bar High School, Class of 2019
USC Viterbi Department of Biomedical Engineering, SHINE 2018



#### Introduction

#### **Atherosclerosis**

- Cardiovascular disease causes over 19 million deaths/year [1].
- Atherosclerotic plaques
  - Build-up of hydroxyapatite (HA), or calcium, fat, and cholesterol [2]
  - Limits blood flow and oxygen supply from the heart to the body





Figure 1. Healthy vs. atherosclerotic artery.
(Source: NIH)

## Objective & Impact of Mentor's Research

Figure 2. Image of a self-assembling micelle (left) and peptide amphiphile (right).

(Source:
ScienceStruck)




My mentor synthesizes nanoparticles called micelles that contain a peptide (HABP) that targets/binds to HA.

- Aids in early theranostics of inflamed plaques before acute cardiac event
- Fluorescent micelle with highest binding affinity to HA was used as imaging probe to locate HA embedded in plaques [2,3]

#### **Skills Learned**

### Hydroxyapatite Binding Peptide (HABP) Conjugation

Hydrophobic tail was linked to HABP [SVSVGMKPSRPC] via thiol-maleimide reaction.



**Figure 3.** Chemical structure of amphiphile DSPE-PEG-Maleimide. (Source: Avanti Polar Lipids, Inc.)

#### **Micelle Synthesis**

Method in which three amphiphiles:

- DSPE-PEG<sub>2000</sub>-HABP (**T**argeting peptide)
- DSPE-PEG<sub>2000</sub>-Methoxy (**E**mpty)
- DSPE-PEG<sub>2000</sub>-FITC (**F**luorescent) were mixed at different molar ratios (see Fig. 4 below) to fabricate micelles exhibiting different binding affinities.

#### **Dynamic Light Scattering (DLS)**

Machine in which lasers are used to measure data such as size and zeta potential of particles.

	Micelle T:E:F Molar	Average Micelle
Sample #	Ratio	Radius (nm)
1	90:0:10	8.8 ± 1.3
2	60:30:10	$6.5 \pm 0.1$
3	45:45:10	5.9 ± 0.2
4	30:60:10	5.2 ± 0.2
5	0:90:10	4.9 ± 0


**Figure 4.** Average micelle radius found using DLS.

#### **Microplate Reader**

A plate reader is used to measure data such as fluorescence and absorbance.

#### **High Performance Liquid Chromatography**

Technique used to purify, separate, and identify different components in a sample. We purified the conjugated HABP sample.



How This Relates to Your STEM
Coursework

**Chemistry:** Stoichiometric calculations, Molarity, Acids and bases, Chemical Structures, Chromatography, Titrations


**Biology:** Cardiovascular System, Peptides, Polymers, Binding sites

# Results from HA Binding Experiments


The binding affinity of HA (powder) to HA targeting micelles were tested.

 Fluorescence intensities of five micelles were found using the Microplate Reader

Total Micelles — Unbound Micelles = Bound Micelles



**Figure 6.** TEM images of fluorescent micelles of Samples 1-5 (A-E). PC: Deborah Chin



**Figure 7.** Results from Binding Experiment with most comparable results. One sample = average of three replicates. Samples 1 and 5 have the highest binding affinity.

## Micelle Synthesis Stoichiometry for Desired Amphiphile Amount:

Volume of desired solvent (mL) x

$$\frac{1 L}{1,000 mL} \times \frac{100 umol}{1 L} \times \frac{1 mol}{1,000,000 umol} \times$$

$$\frac{Molar\ Mass\ of\ PEG\ (g)}{1\ mol}\ x\ \frac{1,000\ mg}{1\ g}$$

= Mass of amphiphile needed (mg)

## Next Steps for You OR Advice for Future SHINE Students

SHINE has given me so much insight on why I want to pursue Biomedical Engineering and to eventually develop therapeutics for diseases.

#### **Advice for future SHINE students:**

- Read the scholarly articles. It helps you understand your research.
- Ask questions!
- Talk to others! Meet new people who will not only help you now, but in the future.

#### References

[1] Naghavi, M., Libby, P., Falk, E., Casscells, S. W., Litovsky, S., Rumberger, J., ... Willerson, J. T. (2003).[2] Lee, J. S., Morrisett, J. D., Tung, C. (2012).[3] Lee, J. S., Tung, C. (2011).

### Acknowledgements

Thank you to Professor Chung for allowing me this amazing opportunity to work in the Chung lab. Thank you to my mentor Deborah, my lab mates (Jonathan, Kelly, Victor, Clarence, Chris, and Yi), and my lab partner Evan for helping me whenever I needed it. And thank you to Dr. Mills and Dr. Herrold for making this experience educational, but fun.